OXYWAY

From serial number 0850000

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A partir del número de serie 0850000

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Description and Operating Instructions

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


Overview

Special information on the unit itself

The symbol on the scale of contents gage **1** is to remind you that the pressure reducer must never be allowed to come into contact with combustible substances (oil, alcohol, grease etc.), as these may cause explosive reactions with oxygen.

Serial number (SN):

08XXXXX
 Serial number
 Year made

The first two digits of the unit number indicate the year of manufacture. The following five digits are used for consecutive numbering of the individual units.

Instructions for use WM 16633 apply to devices up to serial number 0849999.

Safety information in these instructions

Safety information in these instructions for use is marked as follows:

Warning!

Warns of risk of injury and potential material damage.

Caution!

Warns of material damage and potentially false therapy results.

Note:

Contains useful tips.

Contents

Overview	4	5. Hygienic preparation	19
Special information on the unit itself	5	5.1 Intervals	19
Safety information in these instructions	6	5.2 Cleaning	19
1. Description	8	5.3 Disinfecting	19
1.1 Purpose	8	5.4 Changing patients	20
1.2 User qualification	8	6. Functional Checks	21
1.3 Function	8	6.1 Intervals	21
2. Safety Information	11	6.2 Checking the system for leaks	22
2.1 Safety provisions	11	6.3 Checking for leaks in the safety valve	23
2.2 Side effects	13	7. Troubleshooting	24
3. Assembly	14	8. Maintenance	25
3.1 Connecting oxygen cylinder	14	8.1 Replacing the filter screw	25
3.2 Disconnecting pressure reducer	15	8.2 Replacing sealing ring	26
4. Operation	16	8.3 Disposal	26
4.1 Preparing for use	16	9. Product, Accessories	27
4.2 Setting flow	16	9.1 Explanation of numbers	27
4.3 After operation	17	9.2 Standard product	27
4.4 Storage	18	9.3 Accessories	27
4.5 Calculating contents	18	10. Technical Data	29
		11. Warranty	31
		12. Declaration of conformity	32

1. Description

1.1 Purpose

The pressure reducers are for administering medical oxygen at reduced pressure for various devices from the medical and therapy fields and for administering oxygen to the patient via a system (e.g. nasal cannula, oxygen mask, humidifier). The pressure reducers are intended for connection to high-pressure oxygen cylinders (up to 20,000 kPa = 200 bar cylinder pressure).

The pressure reducers are not intended to drive surgical tools.

1.2 User qualification

Patients must be instructed in the use of the device by qualified personnel before therapy begins.

1.3 Function

During use, the oxygen from the cylinder passes through the cylinder valve and connecting union **2** into the body **3** of the pressure reducer. The cylinder pressure (inlet pressure) can be read on the contents gage **1**.


A filter screw with an integrated sintered filter removes any impurities. Another sintered filter reduces the flow of oxygen. The filter screw **11** is located in the opening of connecting union **2**. The second sintered filter is installed between connecting union **2** and body **3**.

In the reduction stage, cylinder pressure is reduced to the operating pressure (outlet pressure) required in each case of between 0.5 and max. 5.5 bar, depending on type. The principle used here, with crater-type bore and spring-loaded piston, ensures particularly constant pressure, despite variations in flow.

All pressure reducers are protected against excessive pressure by a safety valve **10**.

The following pressure reducers are described in these operating instructions:


Fixed pressure reducer (Fix)

A constant flow is generated by a fixed orifice in oxygen outlet **5**.

The pressure reducers listed above differ in the length and position of the connecting union **2**.


Variable pressure reducer (Fine)

The piston is acted on by a second spring which is adjustable by means of a spindle system. This makes it possible to set different operating pressures (outlet pressures). In conjunction with the fixed orifice in oxygen outlet **5**, you can therefore set different flow rates. The flow setting is read off at flow gage **6**.

The pressure reducers listed above differ in the length and position of the connecting union **2**.

Multistep pressure reducer (Fast)

As in the fixed pressure reducer, these produce a constant operating pressure (outlet pressure).

Multistep adjuster **4** houses an orifice disk with orifices of different diameters which are brought into line with oxygen outlet **5** by means of a locating system. This enables you to set different flow rates.


The pressure reducers listed above differ in the length and position of the connecting union **2**.

Oxygen outlets

Oxyway pressure reducers are available with different oxygen outlets: a pressure outlet and a flow outlet (see also "10. Technical Data" on page 29).

Flow outlets are intended for the direct supply of the patient, e.g. via nasal cannula, mask or intermediate humidifier. Flow outlets have a 9/16-18 UNF thread.

Pressure outlets are intended to supply ventilators such as the Medumat, modules or distributor bars. Pressure outlets have a G 3/8" thread.

2. Safety Information

2.1 Safety provisions

For your own safety, the safety of your patients, and to comply with the requirements of EU Directive 93/42 EEC, please observe the following points:

General

Caution!

- Read these operating instructions carefully. They are part of to the pressure reducer and must be available at all times.
- Before working with the pressure reducer, you should understand how it works. In the event of queries, contact your authorized dealer or the manufacturer.
- Use the pressure reducer for the designated purpose only (see "1.1 Purpose" on page 8).
- Have all inspection and servicing work carried out by the manufacturer (WEINMANN) or by a trained authorized dealer.
- Malfunctions and a lack of biocompatibility may result if third-party articles are used. Please bear in mind that in these cases any warranty entitlement and liability shall lapse where the accessories recommended in the instructions for use or original spare parts are not utilised.
- A spare unit must always be kept handy in case of failure.
- Please observe the section "5. Hygienic preparation" on page 19 in order to avoid infection or bacterial contamination.
- **Inadequate supply of oxygen possible.**
Never use pressure reducers with a flow outlet to operate medical devices. The supply pressure of the flow outlets is inadequate for this and could result in malfunctions or failures of the medical devices and thus in a risk to the patient.
- **Device failure possible.**
Read the technical data before connecting a medical product. The use of incompatible devices can lead to the failure of these devices.
- **Inadequate supply of oxygen possible.**
Ensure that the multistep adjuster "engages" properly in the desired setting on OXYWAY Fast pressure reducers. If the flow adjuster is between

two adjacent settings, no oxygen will flow out. This can lead to a risk to the patient.

- **Inadequate supply of oxygen possible.**

Keep the flow in OXYWAY Fine pressure reducers under regular observation and adjust it again if necessary. If the device is used for a prolonged period without interruption, the set flow rate will gradually fall as cylinder pressure drops. This can lead to a risk to the patient.

- **Uncontrolled oxygen output possible!**

The pressure outlets of the pressure reducers have no nonreturn valve to seal the pressure outlet automatically when the connector is removed. When changing devices, therefore, shut off the oxygen cylinder and vent the system (see "4.3 After operation" on page 17) before removing connected devices or tubes.

Oxygen

Warning!

Highly compressed oxygen in conjunction with combustible substances (grease, oil, alcohol etc.) may give rise to spontaneous explosive reactions:

- Keep the units and all screwed unions absolutely free from oil and grease.
- You must never on any account place the pressure reducer in a disinfectant solution or any other cleaning solution. Liquid must never be allowed to enter the pressure reducer because of the resulting explosion risk. A clean cloth, moistened if necessary with clean warm water, is enough for cleaning.
- Be sure to wash your hands before working on the oxygen supply.

Warning!

- If oxygen is being supplied to the respiratory flow, smoking and naked flames are forbidden. **Risk of fire.** The oxygen can accumulate in clothing, bed linen or hair. It can only be removed by thorough ventilation.
- When assembling the unit and when changing cylinders, tighten all screwed unions on the oxygen cylinder and the pressure reducer by hand only. Never use tools. Overtightening damages the screw threads and seals, resulting in leaks.
- Secure the oxygen cylinders so that they cannot fall over. If a cylinder falls on the pressure reducer or valve, these could break off, causing a violent explosion.

Caution!

- Always open the cylinder valve slowly to prevent pressure hammer in the system. Pressure hammer can cause damage to the pressure reducer and the connected devices, or increase the likelihood of spontaneous ignition.
- Do not empty the cylinder completely, as this may allow moist ambient air to enter the cylinder and cause corrosion.

2.2 Side effects

Ask your physician or pharmacist about any side effects which may occur during the inhalation of medical oxygen and read the “medical oxygen” information leaflet enclosed with the product.

3. Assembly

Warning!

Wash your hands thoroughly before doing any work on the oxygen supply. Hydrocarbon compounds (e.g. oil, grease, cleaning alcohol, hand cream or adhesive plaster) can cause explosive reactions if they come into contact with highly compressed oxygen.

Never use wrenches or other tools to tighten or unscrew the union nut.

Warning!

After assembly you must perform a functional check (see “6. Functional Checks” on page 21), to ensure reliable operation.

3.1 Connecting oxygen cylinder


1. Briefly open the valve of the new oxygen cylinder, then shut it again. This is to blow away any particles of dust.

Warning!

In doing so, it is essential to hold the cylinder so that the connector never points at you or another person, otherwise rapidly-dispersing particles of dust could get into the eye and lead to damage. It is also important to prevent O₂ getting into clothing.

2. Make sure that there is absolutely no dirt on the connecting threads of the oxygen cylinder and the pressure reducer.
3. Check that the filter screw **11** is clean. If necessary, fit a new filter screw (see “8.1 Replacing the filter screw” on page 25).

4. Screw the pressure reducer to the cylinder valve with the corrugated union nut of the connecting union **2**. Tighten the union nut by hand.

As soon as the oxygen cylinder is opened, the pressure inside the cylinder pushes the sealing ring into the wedge-shaped annulus between the oxygen cylinder valve and the pressure reducer, thus providing an effective seal.

5. To find out how to connect the device that you want to supply with oxygen, please consult the operating instructions for the relevant device.

3.2 Disconnecting pressure reducer

1. Close the oxygen cylinder valve.
2. **Release the pressure in the system**, e.g. by switching a connected device on and then off or by opening a shut-off valve.

While the system is under pressure, the sealing ring is pressed into the wedge-shaped circular gap between the oxygen cylinder valve and connecting union **2**. For this reason it is impossible to remove the pressure reducer under pressure.


3. Unscrew the union nut **2** by hand.

4. Operation

4.1 Preparing for use

Caution!

Inadequate supply of oxygen possible. Never use pressure reducers with a flow outlet to operate medical devices. The supply pressure of these pressure reducers is not adequate for this and could result in malfunctions or failures of the medical devices and thus in a risk to the patient.


1. **Slowly** open the valve of the oxygen cylinder. Do not turn the valve more than one turn to the left. The contents gage **1** now indicates the cylinder pressure.
2. If necessary, calculate the amount of oxygen left in the cylinder (see "4.5 Calculating contents" on page 18), to ensure that the connected devices can operate long enough.

4.2 Setting flow

Fixed pressure reducer (Fix)

The flow rate is set and cannot be altered.


Variable pressure reducer (Fine)


Set the desired flow (outlet volume) on adjuster knob **7**. You can read the set value on the flow gage **6**.

Caution!

Inadequate supply of oxygen possible. Keep the flow under regular observation and adjust it again if necessary. If the device is used for a prolonged period without interruption, the set flow rate will gradually fall as cylinder pressure drops. This can lead to a risk to the patient.

Multistep pressure reducer (Fast)

Set the multistep adjuster **4** to the desired value.

Caution!

Inadequate supply of oxygen possible.


Ensure that the multistep adjuster “engages” properly in the desired setting on OXYWAY Fast pressure reducers. If the flow adjuster is between two adjacent settings, no oxygen will flow out. This can lead to a risk to the patient.

4.3 After operation

Caution!

Never empty the oxygen cylinder completely. Always ensure that there is a certain residual pressure in the cylinder when you return it for filling, as this prevents moist ambient air from entering the cylinder and causing corrosion.

1. Check the remaining oxygen content on the contents gage **1**. You should change the cylinder in good time, e.g. when the pressure drops below 50 bar, to ensure a sufficiently long operating time.

2. Close the valve on the oxygen cylinder by turning it to the right.
3. In the case of variable or multistep pressure reducers, set the flow to "0". This will prevent pressure hammer in the fittings when the oxygen cylinder is opened again later.

Caution!

The pressure reducer is not a shut-off valve. Even if it is completely closed, it is still possible for a small quantity of oxygen to escape. You should therefore always be sure to close the oxygen cylinder valve after use.

4.4 Storage

If the pressure is not going to be used for a considerable length of time, we recommend the following procedure:

1. Clean the pressure reducer (see "5. Hygienic preparation" on page 19).
2. Store the pressure reducer in a dry place.

Caution!

Be sure to observe the maintenance intervals for stored units as well, as otherwise the unit must not be used when removed from storage.

4.5 Calculating contents

Oxygen content = cylinder volume x cylinder pressure.

	Cylinder volume	x Cylinder pressure	= Oxygen content
Example 1	10 l	x 200 bar	= 2000 l
Example 2	2 l	x 100 bar	= 200 l

5. Hygienic preparation

5.1 Intervals

Clean the device as described in Section “5.2 Cleaning” as required, e.g. after use in a dusty environment.

Disinfect the device as described in Section “5.3 Disinfecting” as required, in particular when changing patients or in the event of direct contact between the patient and the device.

5.2 Cleaning

Warning!

Risk of fire. Do not use cleaning agents. Cleaning agents containing alcohol or moisturizers form a flammable mixture in combination with oxygen.

Caution!

Device damage possible. Ensure that no liquids get into the device, otherwise there may be damage to the device.

- If external cleaning is required, use only a clean cloth. The cloth can be dry or moistened with clean water.

5.3 Disinfecting

Warning!

Risk of injury. Ensure that no disinfectant gets into the tubes, otherwise the disinfectant will get into the patient’s respiratory tract together with the oxygen.

Device damage possible. Ensure that no liquids get into the device, otherwise there may be damage to the device.

- For disinfecting by wiping, we recommend a 0.5 % solution of TERRALIN. You should also follow the instructions for use enclosed with the disinfectant.

5.4 Changing patients

- Perform a hygiene treatment as described in “5.2 Cleaning” and “5.3 Disinfecting”.
- Replace disposable items such as nasal cannula, spiral tube or humidifier with new parts.

6. Functional Checks

After every change of cylinder it is essential to perform a functional check on the pressure reducer.

If this functional check reveals any faults, you must not use the pressure reducer.

A full functional check comprises:

- Visual inspection for mechanical damage;
- "6.2 Checking the system for leaks" on page 22;
- "6.3 Checking for leaks in the safety valve" on page 23.

We recommend you always have the following spare parts available:

- Replacement seal set WM 1148;
- Filter screw WM 30905

6.1 Intervals

After every change of cylinder, and at least every 6 months if not used in the meantime


- Perform a functional check.

Immediately if dirty, at least every two years

- Change filter screw **11** in connecting union **2** (see "8.1 Replacing the filter screw" on page 25).

6.2 Checking the system for leaks

1. Close the system **behind** the pressure reducer, e.g. by switching off all connected devices or by closing the valves on the distribution unit.
2. Check that all screwed unions and tube connections are firmly attached. If necessary, tighten them by hand.

3. **Slowly** open the valve of the oxygen cylinder. You can now read the cylinder pressure on the contents gage **1** of the pressure reducer. For example, a reading of 200 bar means that the cylinder is full, a reading of 100 bar that it is half full.

You should change the cylinder in good time, e.g. when the pressure falls below 50 bar, to ensure a sufficiently long operating time.

4. Close the cylinder valve again.

5. Observe the needle of the contents gage **1** for about 1 minute. If the position of the needle stays constant, the system is free from leaks. If the needle falls steadily, there is a leak in the system.

Rectifying leaks

We recommend you keep a stock of replacement seals for the connectors.

1. Prepare a soapy solution using unperfumed soap.
2. Open the oxygen cylinder.
3. Wet all the screwed unions with the solution. If bubbles form, this indicates a leak.
4. Release the pressure in the system:

Close the oxygen cylinder valve. Briefly switch on the connected device until the contents gage **1** reads "0". Then switch off the connected device.


Caution!

The screwed unions of the oxygen pipes must be tightened only by hand.

5. If there is a leak, remove the faulty seals and replace them with new ones (see "8.2 Replacing sealing ring" on page 26).
6. Then check for leaks again.
7. If the leak cannot be rectified, the unit must be repaired.

6.3 Checking for leaks in the safety valve

1. **Slowly** open the valve on the oxygen cylinder.
2. Set the variable pressure reducers to the maximum flow.

3. Check with your finger whether oxygen is escaping at the safety valve **10**.

If oxygen is escaping, you must have the pressure reducer repaired by the manufacturer (WEINMANN) or by a trained and expert authorized dealer.

7. Troubleshooting

Fault	Cause	Remedy
Leak at cylinder connection	Damaged sealing ring	Replace sealing ring (8.2, page 26)
Leak at device connection	Damaged sealing ring	Replace sealing ring (8.2, page 26)
Leak at safety valve	Excessive pressure rise in pressure reducer	Have unit repaired by manufacturer or trained authorized dealer
Mechanical damage, e.g. to gage or sealing surface		

8. Maintenance

Maintenance work such as inspections or repairs must be carried out by the manufacturer (WEINMANN) or by a trained and authorized dealer.

You must have the pressure reducer serviced regularly as follows.

Immediately if dirty, at least every two years:

- Replace the filter screw **11** (see “8.1 Replacing the filter screw” on page 25).

Every 4 years:

- Replace all parts subject to wear and tear (filter screw, sintered filter, seal elements).
- Perform a final check in accordance with the test instructions.

Oxygen cylinders


Oxygen cylinders must be subjected to retesting at regular intervals by the responsible testing organization. Before using an oxygen cylinder check that the date for retesting has not expired. Oxygen cylinders must not be used if this date has expired. See the data stamped on the cylinder for the date of retesting.

8.1 Replacing the filter screw

You can do this yourself if you are technically minded, or have it done by your authorized dealer.

You must change the filter screw every two years, or if it is clogged.


1. Be sure to wash your hands before changing the filter screw.
2. To unscrew the filter screw you will need a 4 mm Allen key (internal hexagon wrench).
3. Clean the hexagon socket of the filter screw with a clean screwdriver.

4. Unscrew the filter screw to the left, but do not screw it right out yet.
5. To prevent any dirt from falling into the pressure reducer, hold it so that the connector faces down. In this position, unscrew the filter screw completely.
6. Clean the Allen key carefully with washing-up liquid, rinse it well with plenty of clean water and dry it thoroughly.
7. Insert the new filter screw WM 30905 with the aid of the Allen key. Lightly tighten the filter screw by turning to the right.

8.2 Replacing sealing ring

1. Be sure to wash your hands thoroughly before changing the sealing ring.

2. Take a small screwdriver and clean it thoroughly with washing-up liquid. Then rinse it with plenty of clean water and dry the screwdriver thoroughly.

The screwdriver must be clean and free from oil and grease to rule out any risk of explosion.

3. Carefully lift out the sealing ring with the cleaned screwdriver.

Caution!

Take care not to damage the groove in which the sealing ring is seated.

4. If necessary, clean the groove for the sealing ring by wiping with a dry cloth. You may also moisten the cloth slightly with a little clean water.
5. Carefully press the new sealing ring into the groove in the connecting union. Never use a tool for this purpose.

8.3 Disposal

Ensure proper disposal of the unit by consulting an authorized and certified scrap recovery firm.

You can find out their address from your environmental officer or your municipal authorities.

9. Product, Accessories

9.1 Explanation of numbers

Order number:

WM 308xx
└───┬─── Design
 └─── Type

The first three digits of the order number represent the construction principle of the pressure reducer, e.g. "308" for pressure reducers adjustable by latching (Fast). The two last digits represent the specific design, e.g. connection variants.

Device number:

08XXXXX
└───┬─── Serial number
 └─── Year made

The first two digits of the device number indicate the year it was made. The following five digits are serial numbers for the individual devices within a type. The unit number is embossed on the rear of the pressure reducer.

9.2 Standard product

- | | |
|-----------------------------------|----------|
| 1. OXYWAY pressure reducer, 1 pc | |
| 2. Safety tag, 1 pc | WM 16211 |
| 3. Operating instructions, 1 pc | WM 66921 |
| 4. Replacement seal set, 1 pc | WM 1148 |
| 5. Replacement filter screw, 1 pc | WM 30905 |

9.3 Accessories

- | | |
|--|----------|
| 1. Nasal cannula, double-lumen, with earhooks, with tube | WM 1929 |
| – Set I: 20 nasal cannulas (20 x WM 1929) | WM 15220 |
| – Set II: 40 nasal cannulas (40 x WM 1929) | WM 15225 |
| – Set III: 100 nasal cannulas (100 x WM 1929) | WM 15221 |
| 2. Nasal cannula, double-lumen, with tube | WM 1925 |
| – Set I: 20 nasal cannulas (20 x WM 1925) | WM 15251 |
| – Set II: 40 nasal cannulas (40 x WM 1925) | WM 15252 |
| – Set III: 100 nasal cannulas (100 x WM 1925) | WM 15222 |
| 3. Nasal cannula, double-lumen, for children, with tube | WM 1921 |
| – Set I: 20 nasal cannulas (20 x WM 1921) | WM 15344 |
| – Set II: 40 nasal cannulas (40 x WM 1921) | WM 15345 |

	– Set III: 100 nasal cannulas (100 x WM 1921)	WM 15346
4.	Nasal cannula, double-lumen, for babies, with tube	WM 1920
	– Set I: 20 nasal cannulas (20 x WM 1920)	WM 15341
	– Set II: 40 nasal cannulas (40 x WM 1920)	WM 15342
	– Set III: 100 nasal cannulas (100 x WM 1920)	WM 15343

10. Technical Data

Name	Fix I, Lateral outlet Fix I	Fix III Fix III left	Fine I	Fine II	Fine III	Fast I	Fast II	Fast III
Dimensions (WxHxD) in mm	70-82x82-94x69	70-90x94x69	119-239x94x100	119-239x94x100	119-199x94x100	70x94x112	70x94x112	70x94x112
High-pressure hand-screwed connection	country-specific thread							
Weight in g	530	530	730	730	730	710	750	710
Initial pressure p_1 at 15 °C	200 - 10							
Nominal outlet pressure p_2 , flow outlet in bar Outlet thread: 9/16-18 UNF	max. 5,5	-	0,5 - 4,5			max. 5,5		
Nominal outlet pressure p_2 , pressure outlet in bar Outlet thread: G 3/8"	-	4,5	-			-	4,5	-
Flow rate, flow outlet in l/min (at 15° C / 1013 hPa)	4±0,8	-	3-15	1,5-8	0,5-3	1±0,3 2±0,4 3±0,45 4±0,6 5±1 6±1,2 9±1,8 12±2,4 15±3	1±0,3 2±0,4 3±0,45 4±0,6 5±1 6±1,2 9±1,8 12±2,4 15±3	0,5±0,15 1±0,3 1,5±0,45 2±0,4 2,5±0,5 3±0,6 4±0,8 5±1 6±1,2
Nominal pressure outlet flow rate in l/min (at a feed pressure of 100 bar)	-	appr. 120	-			-	90	-
Flow rate, pressure outlet in l/min (at 15° C / 1013 hPa and $p_1 = 20-200$ bar)	-	min. 110	-			-	min. 90	-
Safety valve set pressure in bar	7,4 bar							

Name	Fix I, Lateral outlet Fix I	Fix III Fix III left	Fine I	Fine II	Fine III	Fast I	Fast II	Fast III
Temperature: – Operation – Storage	-20 °C to +60 °C -20 °C to +70 °C							
Product category according to 93/42/EEC	IIb							
Standards satisfied	DIN EN ISO 10524-1							


11. Warranty

- WEINMANN offers a warranty that the product, when used in accordance with requirements, will remain free from defects for a period of two years from date of purchase. For products whose durability is clearly indicated as less than two years, the warranty expires on the expiration date indicated on the packaging or in the user's manual.
- Claims against the warranty can be made only when accompanied by the sales receipt, which must show salesperson and date of purchase.
- We offer no warranty in the case of:
 - Disregard of usage instructions
 - Operating errors
 - Improper use or improper handling
 - Third-party intervention by non-authorized persons for the purpose of device repair
 - Acts of God, e.g., lightning strikes, etc.
 - Transport damage as a result of improper packaging of returned items
 - Lack of maintenance
 - Operational and normal wear and tear, which includes, for example, the following components.
 - Filters
 - Batteries / rechargeable batteries
 - Articles for one-time usage, etc.
 - Failure to use original spare parts.
- WEINMANN is not liable for consequential harm caused by a defect if it is not based on intention or gross negligence. WEINMANN is also not liable for minor physical injury to life or limb resulting from negligence.
- WEINMANN reserves the right to decide whether to eliminate defects, to deliver a defect-free item or to reduce the purchase price by a reasonable amount.
- If WEINMANN rejects a claim against the warranty, it assumes no expense for transport between customer and manufacturer.
- Implied warranty claims remain unaffected by these changes.

12. Declaration of conformity

WEINMANN Geräte für Medizin GmbH + Co. KG declares herewith that the product complies fully with the respective regulations of the Medical Device Directive 93/42/EEC. The unabridged text of the Declaration of Conformity can be found on our website at www.weinmann.de

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Geräte für Medizin GmbH+Co.KG