



## HEINE HSL 150® Hand-held Slit Lamp

### For the examination of the anterior segment



Compact and lightweight, the portable HSL 150 is ideal for situations where a full-size slit lamp is not available or convenient. For the examination of the anterior segment in Human- and Veterinary medicine.

- ∴ **Slit size adjustable** from 0.2mm x 10mm to 4 mm x 14mm.
- ∴ **6x magnification.**
- ∴ **Blue interference filter (FITC)**, can be selected for corneal examination.
- ∴ **Multi-coated optics for maximum light transmission.**
- ∴ **2.5V or 3.5V XHL Xenon Halogen Technology** for bright, white light comparable with the brightness of a classic slit lamp.
- ∴ **The HSL 150 instrument head weighs just 70g**, the lightest slit lamp of its kind.

HSL 150 for BETA or BETA SLIM handles	2.5V	3.5V
HSL 150, head only	<b>C-001.14.602</b>	<b>C-002.14.602</b>
XHL Xenon Halogen spare bulb	<b>X-001.88.098</b>	<b>X-002.88.099</b>
Replacement filter attachment for HSL 150		<b>C-000.14.605</b>

∴ Continuous slit adjustment

∴ Blue filter

∴ Ideal for domiciliary visits

## HEINE HSL 10x® Loupe Attachment

### Loupe attachment with 10x magnification



For HSL 150. With precision magnification optics for crisp, clear images.

- ∴ **10x magnification** when combined with HSL 150 Loupe.
- ∴ **Increased working distance** from examiner to patient.
- ∴ **Large field of view.**
- ∴ **With an additional detachable TP-elastomer Eyepiece** to steady the instrument on the brow and avoid stray light. Can be replaced with a soft TP-elastomer cap, suitable for spectacle-wearers.

HSL 10x Loupe Attachment	<b>C-000.14.606</b>
--------------------------	---------------------

## HEINE HSL 150® Hand-held Slit Lamp Sets

Set complete with HSL 150 and **BETA SLIM handle** (compact) with:



Hand-held Slit Lamp Sets	2.5V	3.5V
BETA SLIM battery handle, in soft pouch, without spare bulb	<b>C-252.10.105</b>	
same, with HSL 10x Loupe Attachment, spare bulb, in hard case	<b>C-268.10.105</b>	
BETA SLIM rechargeable handle and NT 300 charger, spare bulb, in hard case		<b>C-266.20.470</b>
same, with HSL 10x Loupe Attachment, in hard case		<b>C-268.20.470</b>

Set complete with HSL 150, spare bulb and **BETA handle** (standard size), in hard case with:



Hand-held Slit Lamp Sets	2.5V	3.5V (NiMH)	3.5V
BETA battery handle	<b>C-265.10.118</b>		
same, with HSL 10x Loupe Attachment	<b>C-267.10.118</b>		
BETA R rechargeable handle		<b>C-265.20.376</b>	<b>C-265.27.376</b>
same, with HSL 10x Loupe Attachment		<b>C-267.20.376</b>	<b>C-267.27.376</b>
BETA TR rechargeable handle and transformer		<b>C-265.20.384</b>	<b>C-265.27.384</b>
same, with HSL 10x Loupe Attachment		<b>C-267.20.384</b>	<b>C-267.27.384</b>
BETA NT rechargeable handle and NT300 charger		<b>C-265.20.420</b>	
same, with HSL 10x Loupe Attachment		<b>C-267.20.420</b>	
BETA L rechargeable handle and NT 300 charger			<b>C-265.29.420</b>
same, with HSL 10x Loupe Attachment			<b>C-267.29.420</b>

: Rechargeable handle with charge status display in the bottom insert.

230V is the standard mains voltage. If a different voltage is required, e. g. 120V, please specify when ordering.