



Terumo INSULIN SYRINGES

For Single use



**The thinnest, most
comfortable answer for
daily use**

Terumo INSULIN SYRINGES

Comfort first !

In a continuing effort to make the diabetic's life more comfortable, Terumo provides the highest quality in insulin syringes. With our 29G needle we want to minimize trauma and discomfort, which is especially important when the patient needs injections on a daily basis.



Non-Dead-Space design (Myjector)

The permanently attached needle limits dead space to as little as 2µl. This results in less waste of costly insulin and facilitates air bubble removal.

LANCET POINT to minimize discomfort

The needle point is carefully designed to reduce penetration resistance and to guarantee smooth, easy insertion with the least discomfort for the patient.

TERUMO Reference	Colour code		Diameter		Needle length	MI	Units	Packaging		Type
			G	Mm				Unit box	Carton	
BS-01H2516			25	0.5	16 mm	1		100	1800	Slip Tip
BS-01H2613			26	0.45	12 mm	1		100	1800	Slip Tip
BS-01N2516			25	0.5	16 mm	1		100	1800	Slip Tip
BS-01N2613			26	0.45	12 mm	1		100	1800	Slip Tip
BS-N1H2713	■	Orange	27	0.4	12 mm	1	100	100	1800	Myjector
BS-N1H2913	■	Orange	29	0.3	12 mm	1	100	100	1800	Myjector
BS-N1N2713	■	Red	27	0.4	12 mm	1	40	100	1800	Myjector
BS-N1N2913	■	Red	29	0.3	12 mm	1	40	100	1800	Myjector
BS=05M2713	■	Orange	27	0.4	12 mm	0.5	100	100	500	Myjector
BS=05M2913	■	Orange	29	0.3	12 mm	0.5	100	100	500	Myjector
BS=30M2913	■	Orange	29	0.3	12 mm	0.3	100	100	500	Myjector


- : sterilized by EtO

= : sterilized by Gamma



Every day, as healthcare professionals, you are working to ensure life and to contribute to the well-being of your patients.

Every day, Terumo Europe provides you with the clinical systems, means and solutions to keep life flowing.



We keep life flowing

TERUMO

TERUMO CORPORATION
44-1, 2-chome
Hatakeya
Shibuya-ku
Tokyo 151-0072
JAPAN
Tel: +81/3 3374.8111
Fax: +81/3 3374.8195

TERUMO EUROPE N.V.
Researchpark Zone 2 Heaardse
Interfleurwaaier 40
B-3000 Leuven
BELGIUM
Tel: +32/16.39.12.11
Fax: +32/16.40.02.49

TERUMO EUROPE N.V.
BENELUX SALES DIVISION
Interfleurwaaier 40
B-3000 Leuven
BELGIUM
Tel: +32/16.39.25.80
Fax: +32/16.39.25.99
THE NETHERLANDS
Tel: +0000/032.00.96
Fax: +0000/032.04.14

LABORATOIRES TERUMO
FRANCE S.A.
1, Place Aiane
Batiment Uranus
Rue Hélène Boucher
75204 Gargancourt Cedex
FRANCE
Tel: +33/1 30 66 13 00
Fax: +33/1 30 43 60 85

TERUMO DEUTSCHLAND GmbH
Hauptstrasse 57
66760 Eschborn
GERMANY
Tel: +49/61 56.00.220
Fax: +49/61 56.00.22.200

TERUMO EUROPE ESPAÑA SL
Avenida de Burgos 16-D
Edificio EUROCOR
Madrid 28005
ESPAÑA
Tel: +34/91 383.00.60
+34/91 383.91.61
Fax: +34/1 767.04.60

TERUMO ITALIA SRL
Via Simona Martini 143/145
I-00147 Roma
ITALIA
Tel: +39/06.510.61.420
Fax: +39/06.503.04.07

TERUMO SWEDEN AB
Nya Varvet
Buggnads 90
SE-425 71 Västra Frölunda
SWEDEN
Tel: +46/31 748.58.80
Fax: +46/31 748.58.80

TERUMO UK Ltd.
3 Unity Grove
Off School Lane
Kingsley Industrial Park South
Kingsley
Mansfield
L24 9GT
UNITED KINGDOM
Tel: +44/151 452 4800
Fax: +44/151 452 4811

© Registered Trademark
Published by Terumo Europe N.V.

www.terumo.com
www.terumo-europe.com